



Quality: Consistency and Control

The integrity of your work is dependent upon controlled experimental protocols. Our responsibility is to provide nutrition, through our Constant Nutrition® program, that enables you to isolate, distinguish, and define breakthroughs. The entire manufacturing process, from ingredient selection to date-coding of packaged products, is designed to deliver the most consistent nutritional diets available. Ingredients are purchased through approved suppliers who understand and adhere to the demands of research diets.

Computerized formulation and near-infrared spectroscopy allows PMI Nutrition International to maintain the most stringent feed ingredient standards in the industry. Ingredients are analyzed daily for nutritional content as well as non-nutritional contaminants. LabDiet® products are then monitored before approval and release to ensure that quality standards are met. As a result, Certified LabDiet® products provide a managed formulation program, Constant Nutrition® program, that enables you to conduct accurate research by minimizing nutrient variability that can sometimes be associated with fixed formulation diets.

Average Protein Content (%) of Dehulled Soybean Meal 12 **LabDiet* 8 47.5 Percent protein in dehulled soybean meal used in LabDiet* products



Contents

Introduction	⊥
Constant Nutrition® You Can Count On	2
Materials and Methods	3
Results	
Protein/Amino Acids	4
Fat	6
Minerals — Calcium/Phosphorus	7
Vitamins	9
Contaminant Analysis	.10
Summary	.12
References Inside Back Co	wer

Introduction

Animal diets used in laboratory studies must provide consistent and proper nutrition. To assist in experimental control, feed manufacturers must ensure that all batches of feed meet approved nutritional profiles. Certified diets for laboratory animals are commonly used in toxicological studies where control of contaminants, such as pesticides and heavy metals, is imperative. Certified LabDiet® products manufactured by PMI are analyzed prior to sale and guaranteed not to exceed established maximum concentrations of key contaminants. PMI controls variables that could adversely affect study animals and possibly lead to inaccurate conclusions. Certified LabDiet® products meet the recommendations of the FDA's Good Laboratory Practices Program.

This report is an update to our 2004 published bulletin "Nutrient and Contaminant Levels in LabDiet® Certified Laboratory Animal Feeds: 10-Year Technical Bulletin" which reports the results of our ongoing efforts to provide diets that are consistent from batch to batch as well as year to year. The data reported in this bulletin includes that from 1993 to 2002 and now 2003 to 2007.



Constant Nutrition® You Can Count On

PMI's manufacturing process for LabDiet® products ensures consistent levels of nutrients and guarantees minimal concentrations of various contaminants. The Constant Nutrition® program, a managed formulation program exclusive to LabDiet®, ensures a consistent level of nutrients by accounting for biological variation of natural feedstuffs verses fixed formulations, which does not account for that inherent variation resulting in variable nutrients from batch to batch. Therefore, managed formulation requires ingredients be assayed daily and formulations refined to ensure constant nutrient delivery. Alternatively fixed formulation does not consider nutrient variation, thus assays may not be conducted to account for nutrient variability.

Manufactured for a wide variety of animals, key nutrients of LabDiet® products such as protein, fat, crude fiber, calcium, phosphorus, and vitamin C (in monkey and guinea pig rations only) are maintained within established ranges. Concentrations of 30 contaminants are maintained below established levels. Certified LabDiet® products are a key component of animal research studies. As a result, it is essential we control unwanted nutrient- and contaminant-induced variables.

Large discoveries are often built upon data generated over years of research. LabDiet® products have played an important role in building solid baseline data for many research areas.





Combined data from 1993 to 2007 demonstrate that the nutritional content and contaminant levels are well controlled in Certified LabDiet® products manufactured using the Constant Nutrition® program by PMI. Fixed formulation is most appropriately reserved for the equine and pet food industries. It can result in unknown, and sometimes radically variable, nutrient levels as demonstrated in the adjacent table.

PMI's Constant Nutrition® program, driven by our managed formulation philosophy, helps to deliver a consistent nutritional package so that researchers can succeed with reaching their goals.

Based on a database maintained by PMI®, the following example shows how widely protein, alone, can vary in a fixed formula:

Fixed-Formula Diet		The ingredients could conceivably contain the following protein levels (%):		
	%		Low	High
Corn	30.0	Corn	6.5	10.0
Wheat	17.0	Wheat	9.0	14.0
Soybean Meal	15.0	Soybean Meal	45.0	50.0
Midds	6.0	Midds	14.5	17.0
Fish	4.0	Fish	60.0	60.0
Oats	4.0	Oats	12.0	12.0
Alfalfa	3.0	Alfalfa	17.0	17.0
Misc	21.0	Misc	-	-
100.0		Calculated protein level in	formula above:	
			14.5	17.3

Materials and Methods

In order to measure nutritional consistency and derivation of limits of our diets, periodic samples of the following Certified LabDiet® products have been taken over the past 15 years: Certified Rodent Diet 5002, Certified Canine Diet 5007, Certified Guinea Pig Diet 5026, Certified Primate Diet 5048, and Certified Rabbit Diet 5322.

Representative feed samples were sent to an independent laboratory for analysis.* Samples were assayed for proximate constants, essential amino acids, vitamins, and macro and trace minerals. The analytical procedures employed were specified by the Association of Analytical Chemists or were procedures that provide equivalent results. The average analyzed values for nutrients were compared with calculated formulation values. The standard deviation of the analysis for each nutrient was tabulated.



^{*}Nestle-Purina Laboratories, 801 Chouteau Avenue, St. Louis, Missouri 63102.

Results

Protein/Amino Acids

The calculated and analyzed crude protein concentrations of Certified Diets from our 10-year bulletin as well as those collected in the past 5 years are provided in Table 1A. Analyzed concentrations of crude protein tend to be slightly greater than that of the calculated values. The consistency of crude protein in certified diets is evident based on the data collected for the past 15 years. Chart 1B indicates the consistency of protein in 5002 during the same period. Ingredients are analyzed daily for crude protein using NIR to ensure our formulations and finished products never fall below the level guaranteed by the tag.

Proteins are composed of amino acids. Laboratory animals have specific dietary requirements for numerous amino acids such as arginine, histidine, isoleucine, leucine, lysine, methionine, phenylalanine, threonine, tryptophan and valine. All of these must be in sufficient concentrations in the diet to ensure proper animal performance (growth, gestation and lactation). Because lysine and methionine are rate limiting amino acids, their concentrations over time are exhibited in Tables 1C and 1D. Similar to the first 10 years, data collected in the past 5 years illustrates a very consistent product from batch to batch.



Table 1A.
Protein Levels in
Certified LabDiet
Products Over
15-year Period

10- and 15-Year Average Protein Levels (%) in Certified LabDiets®						
DIET	Calculated Value %	# Samples 10 Years	10-Year Average	# Samples 15 Years	15-Year Average	2003-2007 3rd 5-year Average
Certified Rodent Diet 5002	20.10	30	21.29+/-0.41	47	21.23+/-0.40	21.09+/-0.35
Certified Canine Diet 5007	25.50	30	26.82+/-0.57	47	26.70+/-0.57	26.39+/-0.45
Certified Guinea Pig Diet 5026	18.50	29	20.13+/-0.47	47	20.07+/-0.51	19.80+/-0.59
Certified Primate Diet 5048	25.00	29	26.96+/-0.61	47	26.81+/-0.61	26.55+/-0.51
Certified Rabbit Diet 5322	16.20	29	17.48+/-0.68	47	17.36+/-0.64	17.08+/-0.41
						+/- Standard Deviation

Chart 1B.
Protein Levels in
Certified Rodent
Diet 5002 Over
15-Year Period






Table 1C.
Lysine Levels in
Certified LabDiet
Products Over a
15-year Period

Lysine % in Certified LabDiets®						
DIET	Calculated Value %	# Samples 10 Years	10-Year Average	# Samples 15 Years	15-Year Average	2003-2007 3rd 5-year Average
Certified Rodent Diet 5002	1.20	30	1.15+/-0.06	47	1.16+/-0.05	1.16+/-0.04
Certified Canine Diet 5007	1.10	30	1.24+/-0.07	47	1.25+/-0.06	1.26+/-0.06
Certified Guinea Pig Diet 5026	1.00	29	1.03+/-0.13	47	1.04+/-0.11	1.04+/-0.06
Certified Primate Diet 5048	1.20	29	1.21+/-0.56	47	1.22+/-0.05	1.22+/-0.04
Certified Rabbit Diet 5322	0.80	29	0.88+/-0.06	47	0.87+/-0.05	0.86+/-0.04
						+/- Standard Deviation

Table 1D.

Methionine Levels in
Certified LabDiet
Products Over a
15-year Period

	Methionine % in Certified LabDiets®							
DIET	Calculated Value %	# Samples 10 Years	10-Year Average	# Samples 15 Years	15-Year Average	2003-2007 3rd 5-year Average		
Certified Rodent Diet 5002	0.40	30	0.41+/-0.02	47	0.42+/-0.02	0.43+/-0.02		
Certified Canine Diet 5007	0.40	30	0.41+/-0.03	47	0.42+/-0.03	0.43+/-0.03		
Certified Guinea Pig Diet 5026	0.40	29	0.39+/-0.07	47	0.38+/-0.05	0.38+/-0.02		
Certified Primate Diet 5048	0.50	29	0.47+/-0.02	47	0.47+/-0.03	0.49+/-0.03		
Certified Rabbit Diet 5322	0.40	29	0.34+/-0.02	47	0.34+/-0.02	0.35+/-0.03		
						+/- Standard Deviation		



Fat

Fat was analyzed using 2 different procedures: ether extract and acid hydrolysis. Ether extraction of fat measures neutral fat and is an appropriate assay for most feeds. However, fat in extruded animal foods, such as canine and primate diets, must be hydrolyzed with dilute acid (acid hydrolysis) to free the fat molecules, which is then extracted and measured. Acid hydrolysis of fat also removes bound fat from ingredients, along with other forms of lipids (i.e. phospholipids). Consequently, fat measured by acid hydrolysis is generally 1 or 2 percentage units greater than fat measured by ether extract.

Calculated and analyzed crude fat levels for LabDiet® 5002 and the extruded primate diet, LabDiet® 5048, are provided in Table 2. Over time, crude fat concentrations have remained relatively constant and tend to be similar to calculated values. Consistent dietary fat is important to maintain steady energy levels and feed consumption.

You may note that calculated fat (ether extract) values do not account for fat bound by the extrusion process. Thus, the analyzed values, by ether extract, for extruded products (Canine, Feline, and Primate) will be less than the calculated values. Although one extraction process works better for various diet forms (pelleted vs. extruded), we continue to use both extraction methods when analyzing for crude fat regardless of form.

Table 2.
Fat Levels in Certified
Rodent Diet 5002
and Certified Primate
Diet 5048 Products
Over 15-Year Period

	Certified Rod	ent Diet 5002 Fat A.H., % ²	Certified Primate Diet 5048 Fat E.E., %1 Fat A.H., %2
10-Year Average	4.61+/-0.45	5.78+/-0.25	3.45+/-0.71 6.35+/-0.62
2003	4.67+/-0.35	5.71+/-0.25	3.69+/-0.20 6.52+/-0.44
2004	4.74+/-0.10	5.61+/-0.29	3.06+/-0.10 6.19+/-0.48
2005	4.46+/-0.10	5.60+/-0.26	4.04+/-0.10 6.23+/-0.80
2006	4.75+/-0.05	5.54+/-0.23	2.86+/-0.21 6.01+/-0.74
2007	5.08+/-0.34	5.85+/-0.27	3.62+/-0.28 5.97+/-0.28
15-Year Average	4.65+/-0.40	5.74+/-0.27	3.60+/-1.13 6.29+/-0.63
Calculated Value, %	4.50	5.10	4.70 5.60
¹ Ether Extract ² Acid Hydr	olysis +/- Standard D	eviation	





Minerals - Calcium/Phosphorus

Calcium and phosphorus are the primary macro minerals in mammals and are necessary for maintenance of bone, teeth, nerve function, energy, metabolism and milk production. Minerals are also cofactors for various enzymes. Calcium and phosphorus must be delivered to laboratory animals at proper dietary concentrations with a calcium:phosphorus ratio maintained at 1:1 or greater.

Observed mineral concentrations are provided in Charts 3A and 3B for diets 5002, 5007 and 5026. Analyzed values for calcium and phosphorous for the past 5 years are consistent with values tabulated and published in the 10-Year Technical Bulletin further demonstrating nutrient consistency due to Constant Nutrition® formulation.

Chart 3A.

Calcium Levels (%) in Certified LabDiet® Products Over 15-Year Period









Chart 3B.
Phosphorus Levels (%) in Certified LabDiet® Products Over 15-Year Period











Minerals - Calcium/Phosphorus (cont.)


10-Yr Avg


Chart 3A.
Calcium Levels (%) in Certified LabDiet® Products Over 15-Year Period


2004

15-Yr Avg

Chart 3B.
Phosphorus Levels(%) in Certified LabDiet® Products Over 15-Year Period











Vitamins

Riboflavin, pantothenic acid, and vitamin B_{12} analyzed in Certified LabDiet® products are provided in Tables 4A, 4B and 4C, respectively. Values tabulated during the first 10 years of our collection period are consistently similar to those collected within the past 5 years. Some variation has been observed when comparing calculated values to analyzed values. Analytical and sampling variation may explain some of the discrepancies seen between these values.

Vitamin Levels in Certified LabDiet Products Over 15-Year Period

Table 4A. Riboflavin

	Riboflavin Levels (PPM)						
DIET	Calculated Value	10-Year Average	15-Year Average				
5002	8.00	7.80+/-1.49	7.81+/-1.30				
5007	4.50	4.61+/-0.84	4.70+/-0.89				
5026	6.00	6.63+/-1.53	6.76+/-1.47				
5048	8.50	8.44+/-1.79	8.55+/-1.62				
5322	5.00	6.44+/-1.81	6.49+/-1.53				
			+/- Standard Deviation				

Table 4B. Pantothenic Acid

Pantothenic Acid Levels (PPM)								
DIET	Calculated Value 10-Year Average 15-Year Average							
5002	17.00	17.19+/-1.97	17.08+/-1.30					
5007	20.00	20.73+/-2.30	20.81+/-3.70					
5026	19.00	21.50+/-1.82	21.61+/-3.27					
5048	50.00	56.34+/-4.67	55.95+/-6.69					
5322	19.00 21.52+/-1.96 22.10+/-2.79							
			+/- Standard Deviation					

Table 4C. Vitamin B₁₂

Vitamin B ₁₂ Levels (MCG/KG)							
DIET	Calculated Value	10-Year Average	15-Year Average				
5002	19.80	26.70+/-8.68	27.34+/-8.32				
5007	26.90	37.57+/-11.95	37.85+/-11.82				
5026	13.20	19.48+/-5.47	20.70+/-6.85				
5048	46.00	48.24+/-8.09	49.57+/-9.29				
5322	6.60	14.83+/-5.13	15.37+/-5.26				
			+/- Standard Deviation				

Contaminant Analysis

Heavy metal and/or chemical contamination of lab animal diets used for GLP studies is unacceptable. Thus, these substances must be closely monitored. Composite samples from each manufacturing lot of Certified LabDiet® product are analyzed for 30 contaminants, including heavy metals, aflatoxin, chlorinated hydrocarbons and organophosphates. The certification profile has a maximum allowable concentration for each contaminant measured. Analyzed concentrations of all 30 contaminants were consistently below maximum allowable levels. Concentrations of arsenic, malathion, lead, and PCBs are exhibited in Tables 5A, 5B, 5C and 5D, respectively. For all contaminant assays, the lower limit of detection

is considerably lower than our guaranteed maximum limits. For many of the contaminants analyzed concentrations fall well below that detection limit. Because of our tight controls, data collected during the 15-year period is consistent and reliable.

In the periodic samples submitted to an independent laboratory for this study, there was no derivation of limits of aflatoxin detected. This indicates that our Certified LabDiet® program for controlling mycotoxin contamination has achieved its intended goal.

Contaminant Analysis of Certified Labdiet® Products Over 15-Year Period

Table 5A.

Arsenic Concentrations
in Certified LabDiet®
Products Over Test
Period of 1993-2007

	Arsenic, ppm						
YEAR	5002	5007	5026	5048	5322		
10-Year Average	0.24+/-0.14	0.27+/-0.16	0.23+/-0.15	0.23+/-0.14	0.21+/-0.15		
2003	0.10+/-0.12	<0.02	0.04+/-0.12	0.07+/-0.18	< 0.02		
2004	0.03+/-0.09	0.04+/-0.13	0.17+/-0.11	0.02+/-0.07	0.15+/-0.14		
2005	0.03+/-0.14	<0.02	0.25+/-0.03	< 0.02	0.23+/-0.17		
2006	<0.02	<0.02	0.19+/-0.14	0.03+/-0.07	0.33+/-0.11		
2007	< 0.02	<0.02	0.29+/-0.16	< 0.02	0.47+/-0.22		
15-Year Average	0.21+/-0.15	0.24+/-0.17	0.22+/-0.15	0.19+/-0.16	0.19+/-0.15		
Max. Conc. Levels	1.00 PPM	1.00 PPM	1.00 PPM	1.00 PPM	1.00 PPM		
Limits of Detection	0.02	0.02	0.02	0.02	0.02		
					+/- Standard Deviation		

Table 5B.

Malathion Concentrations
in Certified LabDiet®
Products Over Test
Period of 1993-2007

Malathion, ppm							
YEAR	5002	5007	5026	5048	5322		
10-Year Average	0.05+/-0.06	0.04+/-0.06	0.04+/-0.05	0.04+/-0.06	0.05+/-0.06		
2003	< 0.02	<0.02	0.02+/-0.02	0.02+/-0.03	0.05+/-0.05		
2004	<0.02	<0.02	<0.02	< 0.02	< 0.02		
2005	< 0.02	<0.02	<0.02	< 0.02	< 0.02		
2006	0.02+/-0.02	<0.02	0.02+/-0.02	< 0.02	0.02+/-0.02		
2007	0.03+/-0.05	<0.02	0.05+/-0.06	< 0.02	0.07+/-0.07		
15-Year Average	0.04+/-0.05	0.03+/-0.05	0.04+/-0.05	0.03+/-0.05	0.04+/-0.06		
Max. Conc. Levels	0.5 PPM						
Limits of Detection	0.02	0.02	0.02	0.02	0.02		
					+/- Standard Deviation		



Contaminant Analysis of Certified Labdiet® Products Over 15-Year Period

Table 5C. Lead Concentrations in Certified LabDiet® Products Over Test Period of 1993-2007

Lead, ppm						
YEAR	5002	5007	5026	5048	5322	
10-Year Average	0.22+/-0.05	0.83+/-0.09	0.27+/-0.08	0.18+/-0.05	0.30+/-0.09	
2003	0.13+/-0.05	0.36+/-0.29	0.18+/-0.07	0.11+/-0.05	0.19+/-0.07	
2004	0.18+/-0.03	0.28+/-0.04	0.26+/-0.04	0.17+/-0.02	0.27+/-0.04	
2005	0.20+/-0.03	0.31+/-0.05	0.25+/-0.05	0.19+/-0.03	0.29+/-0.05	
2006	0.19+/-0.03	0.29+/-0.04	0.28+/-0.03	0.18+/-0.02	0.28+/-0.04	
2007	0.20+/-0.03	0.31+/-0.04	0.32+/-0.05	0.21+/-0.03	0.31+/-0.04	
15-Year Average	0.21+/-0.05	0.68+/-0.25	0.73+/-0.08	0.18+/-0.05	0.29+/-0.08	
Max. Conc. Levels	1.5 PPM					
Limits of Detection	0.05	0.05	0.05	0.05	0.05	
					+/- Standard Deviation	


Table 5D.
PCBs Concentrations
in Certified LabDiet®
Products Over Test
Period of 1993-2007

PCBs, ppm					
YEAR	5002	5007	5026	5048	5322
10-Year Average	< 0.15	<0.15	< 0.15	< 0.15	< 0.15
2003	< 0.15	< 0.15	< 0.15	< 0.15	< 0.15
2004	< 0.15	<0.15	< 0.15	< 0.15	< 0.15
2005	< 0.15	<0.15	< 0.15	< 0.15	< 0.15
2006	< 0.15	<0.15	< 0.15	< 0.15	<0.15
2007	< 0.15	<0.15	< 0.15	< 0.15	< 0.15
15-Year Average	< 0.15	<0.15	< 0.15	< 0.15	<0.15
Max. Conc. Levels	0.15 PPM				
Limits of Detection	0.15	0.15	0.15	0.15	0.15
					+/- Standard Deviation



Summary

The cumulative results of what is now a 15-year collection period demonstrates that nutrient and contaminant concentrations are consistent with our Certified LabDiet® products. LabDiet® has placed strict constraints on our certified products through the use of our managed formulation program of Constant Nutrition®, approved suppliers, daily nutrient analyses, and defined manufacturing protocols. In the end, this translates to greater control and more consistency with your research experiments.





References

- 1. Cohen LA, Epstein M, Sao Pabon V, Meschter C, Zang E. Interactions between 4-HPR and diet in NMU-induced mammary tumorigenesis. *Nutr Cancer*. 1994;21:271-283.
- 2. National Research Council. *Nutrient Requirements of Laboratory Animals*. Washington, DC: National Academy Press; 1995:112-113.
- 3. Verchere CB, D'Alessio DA, Palmiter RD, et al. Islet amyloid formation associated with hyperglycemia in transgenic mice with pancreatic beta cell expression of human islet amyloid polypeptide. *Proc Natl Acad Sci USA*. 1996;93:3492-3496.
- 4. Cockell KA, L'Abbé MR, Belonje B. The concentrations and ratio of dietary calcium and phosphorus influence development of nephrocalcinosis in female rats. *J Nutr.* 2002;132:252-256.
- 5. Ritskes-Hoitinga J, Lemmens AG, Danse LHJC, Beynen AC. Phosphorus-induced nephrocalcinosis and kidney function in female rats. *J Nutr.* 1989;119:1423-1431.
- 6. Stonard MD, Samuels DM, Lock EA. The pathogenesis and effect on renal function of nephrocalcinosis induced by different diets in female rats. *Food Chem Toxicol*. 1984;22:139-146.
- 7. Anand CR, Linkswiler HM. Effect of protein intake on calcium balance of young men given 500 mg calcium daily. *J Nutr.* 1974;104:695-700.



Committed to providing premium products to research professionals around the world, LabDiet® conducts innovative laboratory animal research to provide Constant Nutrition® that will make your research a success.





For more information, visit www.labdiet.com or call 1-800-227-8941. LabDiet® is a registered trademark of PMI Nutrition International LLC.